

## Rat IgG2a, Lambda Isotype Control [B39-4] {PerCP/Cyanine5.5} (STJA0002758)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	PerCP/Cyanine5.5
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG2a, lambda
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

### ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---