

Rat IgG2a, Lambda Isotype Control [B39-4] {PerCP/Cyanine5.5} (STJA0002774)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	PerCP/Cyanine5.5
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG2a, lambda
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---