

## Anti-laminin Alpha 4 antibody [3H2] (STJA0002806)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/IHC/WB/IP
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G purification
Dilution Range	ELISA titer 0.02 $\mu$ g/ml; WB, reduced conditions; 0.5 to 1 $\mu$ g/ml; IP 0.5 mg/ml; IHC dilution 5 to 10 $\mu$ g/ml Monoclonal antibody working titer has to be established practically for each particular antigen and assay format
Formulation	PBS pH 7.4, with 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at -20 to -70 °C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles

### TARGET INFORMATION

Gene ID	<a href="#">3910</a>
Gene Symbol	<a href="#">LAMA4</a>
UniProt ID	<a href="#">LAMA4_HUMAN</a>
Specificity	Human laminin Alpha 4

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

