

Anti-DNase I antibody (STJA0002809)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Rabbit
Reactivity	DNASE1

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Bovine DNaseI-affinity purification
Dilution Range	ELISA 0.05 to 0.1 µg/ml; WB 0.1 to 0.2 µg/ml; IF 0.2 to 12 µg/ml Antibody working titer has to be established practically for each particular antigen and assay format.
Formulation	Concentrated ammonium sulphate in PBS pH 7.4.
Isotype	IgG
Storage Instruction	Store at +4°C upon receipt. As product is (NH ₄) ₂ SO ₄ precipitate, mix well by pipetting or vortexing prior use. Do not freeze

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---