

Anti-Papillomavirus type 1 E2-protein antibody (182-190 aa) [1E2] (STJA0002810)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/IF/IP
Host / Source	Mouse
Reactivity	PV

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A purification
Dilution Range	ELISA 0.2 to 1 μ g/ml; WB 0.5 to 2 μ g/ml; IF 0.33 to 1 μ g/ml Antibody working titer has to be established practically for each particular antigen and assay format
Formulation	PBS pH 7.4, 50% glycerole
Isotype	IgG
Storage Instruction	Store at -20 or -70 °C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles as product degradation may result

TARGET INFORMATION

Immunogen Region	182-190 aa
------------------	------------

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

