

Anti-Papillomavirus type 1 E2-protein antibody (197-208 aa) [3F12] (STJA0002814)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/IF
Host / Source	Mouse
Reactivity	PV

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A purification
Dilution Range	ELISA 0.05 to 0.1 $\mu\text{g/ml}$; WB 0.33 to 1 $\mu\text{g/ml}$; IF 1 $\mu\text{g/ml}$ Monoclonal antibody working titer has to be established practically for each particular antigen and assay format.
Formulation	PBS with 50% glycerol
Isotype	IgG
Storage Instruction	Store at -20 or -70 °C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles as product degradation may result.

TARGET INFORMATION

Immunogen Region	197-208 aa
------------------	------------

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

