

Anti-HCV subtype 1b NS5B protein antibody (77-86 aa) [7F12] (STJA0002828)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/WB/IP
Host / Source	Mouse
Reactivity	Human/HCV

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G purification
Dilution Range	WB in reduced conditions 1 μ g/ml; Optimal conditions for IP should be determined for each particular application. Monoclonal antibody working titer has to be established practically for each particular antigen and assay format.
Formulation	PBS pH 7.4, with 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at -20 to -70 °C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles as product degradation may result.

TARGET INFORMATION

Immunogen Region	77-86 aa
Specificity	Human, HCV subtype 1b NS5B, Epitope mapped to amino acids 77-86 (TVKAKLLSVE)

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---