

Anti-HIV-1 Tat antibody (71-81 aa) [N1] (STJA0002848)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ChIP
Host / Source	Mouse
Reactivity	HIV-1

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G purification
Dilution Range	ELISA 0.08 to 0.4 μ g/ml; WB 2 μ g/ml Monoclonal antibody working titer has to be established practically for each particular antigen and assay format
Formulation	PBS pH 7.4, with 0.1 % sodium azide
Isotype	IgG3
Storage Instruction	Store at -20 to -70 °C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles

TARGET INFORMATION

Immunogen Region	71-81 aa
Specificity	Reacts with HIV-1 Tat of HXB2, BH10, IIB and HAN strains. Mapped to amino acids 71-81 (KQPASQPRGDP – B subtype a.a. consensus)

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

