

Anti-HIV-1 Tat antibody (59-67 aa) [N8] (STJA0002850)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ChIP
Host / Source	Mouse
Reactivity	HIV-1

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G purification
Dilution Range	ELISA 0.05 to 0.1 μ g/ml; WB 2 μ g/ml. Monoclonal antibody working titer has to be established practically for each particular antigen and assay format
Formulation	PBS pH 7.4, with 0.1% sodium azide
Isotype	IgG2b
Storage Instruction	Store at -20 to -70 °C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles

TARGET INFORMATION

Immunogen Region	59-67 aa
Specificity	Reacts in ELISA and WB with HIV-1 Tat of HXB2, BH10, IIIB and HAN strains. Mapped to amino acids 59-67 (PQDSQTHQV – B subtype a.a. consensus)

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

