

Anti-HPV-11 E7 antibody (1-35 aa) [8G6] (STJA0002854)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/WB/IF/IP
Host / Source	Mouse
Reactivity	HPV-11

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G purification
Dilution Range	ELISA 0.045 to 0.09 $\mu\text{g/ml}$ WB 0.9 $\mu\text{g/ml}$ IP 1.8 $\mu\text{g/ml}$ Monoclonal antibody working titer has to be established practically for each particular antigen and assay format
Formulation	PBS pH 7.4, with 0.1% sodium azide
Isotype	IgG2a
Storage Instruction	Store at -20 to -70°C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles

TARGET INFORMATION

Immunogen Region	1-35 aa
Specificity	Reacts with HPV-11 E7-protein. Mapped to amino acids 1-35

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---