

Anti-HPV-11 E7 antibody (36-70 aa) [9H5] (STJA0002855)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	ELISA/WB/IF/IP
	Host / Source	Mouse
	Reactivity	HPV-11
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Conjugation	Unconjugated
	Purification	Protein G purification
	Dilution Range	ELISA 0.05 to 0.1 μ g/ml; WB 1 μ g/ml; IP 2 μ g/ml. Monoclonal antibody working titer has to be established practically for each particular antigen and assay format
	Formulation	PBS pH 7.4, with 0.1% sodium azide
	Isotype	IgG2a
	Storage Instruction	Store at-20 to-70°C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles
TARGET INFORMATION		
	Immunogen Region	36-70 aa
	Specificity	Reacts with HPV-11 E7-protein. Mapped to amino acids 36-70
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.