

Anti-HPV-18 E2 antibody (1-83 aa) [2E7] (STJA0002856)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/WB/IF
Host / Source	Mouse
Reactivity	HPV-18

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G purification
Dilution Range	ELISA 0.1 μ g/ml, WB 1 μ g/ml This antibody does not work in IP. Monoclonal antibody working titer has to be established practically for each particular antigen and assay format
Formulation	PBS pH 7.4, with 0.1% sodium azide
Isotype	IgG
Storage Instruction	Store at -20 to -70°C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles

TARGET INFORMATION

Immunogen Region	1-83 aa
Specificity	Reacts with HPV-18 E2-protein. Mapped to amino acids 1-83

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---