

Anti-MANF antibody [1G12] (STJA0002862)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G purification
Dilution Range	ELISA 0.1 to 0.2 μ g/ml; WB 0.5 to 2 μ g/ml; IF 1 to 10 μ g/ml Monoclonal antibody working titer has to be established practically for each particular antigen and assay format
Formulation	PBS pH 7.4, with 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at -20°C to -70°C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles

TARGET INFORMATION

Gene ID	7873
Gene Symbol	MANF
UniProt ID	MANF_HUMAN
Specificity	Human MANF

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---