

## Anti-Hev b1 antibody [1-1] (STJA0002869)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Mouse
Reactivity	HEVB1

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G purification
Dilution Range	ELISA 132 000. Reacts as detection antibody in Capture ELISA with Hev b1 binding monoclonal antibody A3-701-100. WB 12000 to 14000 Monoclonal antibody working titer has to be established practically for each particular antigen and assay format
Formulation	PBS, 0.1% sodium azide
Isotype	IgG2a
Storage Instruction	Store at -20 to -70 °C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles as product degradation may result

### TARGET INFORMATION

Specificity	Recombinant and native Hev b1
-------------	-------------------------------

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

