

Anti-Hev b5 antibody [4] (STJA0002872)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | ELISA |
| Host / Source | Mouse |
| Reactivity | HEVB5 |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Clonality | Monoclonal |
| Conjugation | Unconjugated |
| Purification | Protein G purification |
| Dilution Range | ELISA 132 000 Reacts as binding antibody in capture ELISA with Hev b5 detection monoclonal antibody A3-708-100. |
| Formulation | PBS, 0.1% sodium azide |
| Isotype | IgG1 |
| Storage Instruction | Store at -20 to -70 °C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles as product degradation may result |

TARGET INFORMATION

| | |
|-------------|-------------------------------|
| Specificity | Recombinant and native Hev b5 |
|-------------|-------------------------------|

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|