

Anti-DNase I antibody [4B10] (STJA0002876)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Mouse
Reactivity	DNASE1

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G purification
Dilution Range	ELISA 0.01 to 0.02 μ g/ml; WB 0.1 to 0.2 μ g/ml (bovine DNase I) Antibody working titer has to be established practically for each particular antigen and assay format
Formulation	PBS pH 7.4, with 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at -20 to -70 °C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles as product degradation may result

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

