

Anti-Neurturin antibody (STJA0002881)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/IF
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Immunogen affinity purification
Dilution Range	WB 0.2 to 1 μ g/ml; ELISA 0.01 to 0.02 μ g/ml; IF 0.5 to 2 μ g/ml Antibody doesn't cross-react to other GDNF-family proteins Polyclonal antibody working titer has to be established practically for each particular antigen and assay format
Formulation	Concentrated ammonium sulphate in PBS pH 7.4
Isotype	IgG
Storage Instruction	Store at +4°C upon receipt. As product is (NH ₄) ₂ SO ₄ precipitate, mix well by pipetting or vortexing prior use

TARGET INFORMATION

Gene ID	4902
Gene Symbol	NRTN
UniProt ID	NRTN_HUMAN
Specificity	Human NRTN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

