

Anti-Neurturin antibody [4G8] (STJA0002883)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G purification
Dilution Range	WB 0.2 to 1 μ g/ml; ELISA 0.05 to 0.1 μ g/ml; IF 2 to 10 μ g/ml Monoclonal antibody working titer has to be established practically for each particular antigen and assay format
Formulation	PBS pH 7.4, with 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at -20 or -70 °C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles as product degradation may result

TARGET INFORMATION

Gene ID	4902
Gene Symbol	NRTN
UniProt ID	NRTN_HUMAN
Specificity	Human NRTN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---