

Anti-GFRa-4 antibody (STJA0002885)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/IF
Host / Source	Chicken
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	GFRa-4 affinity purified
Dilution Range	WB 0.1 to 0.2 μ g/ml; ELISA 3 to 6.7 ng/ml; IF 0.3 to 1 μ g/ml Polyclonal antibody working titer has to be established practically for each particular antigen and assay format
Formulation	PBS pH 7.4
Isotype	IgY
Storage Instruction	Store at -20 to -70°C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles.

TARGET INFORMATION

Gene ID	64096
Gene Symbol	GFRA4
UniProt ID	GFRA4_HUMAN
Specificity	Human GFRa-4, doesn't cross-react with hGFRa-1, hGFRa-2, hGFRa-3

ADDITIONAL INFORMATION

Note: STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

