

Anti-Recombinant-SARS-CoV-2 Nucleocapsid protein antibody [16C11] (STJA0002888)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Human
Reactivity	SARS-CoV-2

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A affinity chromatography following desalting
Dilution Range	ELISA 0.08-0.16 µg/ml
Formulation	PBS pH 7.4
Isotype	IgG1 lambda
Storage Instruction	The product is stable at 4 °C for several days. Store at -20 °C to -70 °C long term. Avoid multiple freeze-thaw cycles.

TARGET INFORMATION

Specificity	SARS-CoV-2 Nucleocapsid protein
-------------	---------------------------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---