

## Anti-Recombinant-SARS-CoV-2 Spike RBD antibody [27B12] (STJA0002892)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA
Host / Source	Human
Reactivity	SARS-CoV-2

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A affinity chromatography following desalting
Dilution Range	Indirect ELISA 0.02 – 0.04 ng/ml
Formulation	PBS pH 7.4
Isotype	IgG1 kappa
Storage Instruction	Store at -20 °C to -70 °C long term. Avoid multiple freeze-thaw cycles.

### TARGET INFORMATION

Specificity	SARS-CoV-2 Spike RBD
-------------	----------------------

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---