

Anti-Artemin antibody [2G7] (STJA0002898)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/IF/IHC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G purification
Dilution Range	ELISA 0.2 to 1 μ g/ml; IF 1 to 20 μ g/ml; IHC 5 μ g/ml. Conformational antibody, not suitable for Western Blot application. Monoclonal antibody working titer has to be established practically for each particular antigen and assay format
Formulation	PBS pH 7.4, with 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at -20 to -70°C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles

TARGET INFORMATION

Gene ID	9048
Gene Symbol	ARTN
UniProt ID	ARTN_HUMAN
Specificity	Human Artemin

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---