

## Anti-Ribonuclease 7 antibody [3C2] (STJA0002899)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | ELISA/WB/IF/IHC    |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |                                                                                                                                                                                                              |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                                                                                                                                   |
| Conjugation         | Unconjugated                                                                                                                                                                                                 |
| Purification        | Protein G purification                                                                                                                                                                                       |
| Dilution Range      | ELISA 0.05 to 0.1 $\mu$ g/ml; WB 0.5 to 1 $\mu$ g/ml; IF 2 to 20 $\mu$ g/ml; IHC 5 to 20 $\mu$ g/ml Monoclonal antibody optimal dilution has to be established practically for each antigen and assay format |
| Formulation         | PBS pH 7.4, with 0.1 % sodium azide                                                                                                                                                                          |
| Isotype             | IgG1                                                                                                                                                                                                         |
| Storage Instruction | Store at -20 °C to -70 °C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles as product degradation may result                                                       |

### TARGET INFORMATION

|             |                              |
|-------------|------------------------------|
| Gene ID     | <a href="#">84659</a>        |
| Gene Symbol | <a href="#">RNASE7</a>       |
| UniProt ID  | <a href="#">RNASE7_HUMAN</a> |
| Specificity | Human Ribonuclease 7         |

### ADDITIONAL INFORMATION

**Note** STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

