

Anti-Ribonuclease 7 antibody [4F9] (STJA0002901)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/IHC/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G purification
Dilution Range	ELISA 0.02 to 0.05 μ g/ml; WB 0.5 to 1 μ g/ml; IF 1 to 10 μ g/ml; IHC 5 to 20 μ g/ml Monoclonal antibody optimal dilution has to be established practically for each antigen and assay format
Formulation	PBS pH 7.4, with 0.1 % sodium azide
Isotype	IgG1
Storage Instruction	Store at -20 °C to -70 °C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles as product degradation may result

TARGET INFORMATION

Gene ID	84659
Gene Symbol	RNASE7
UniProt ID	RNASE7_HUMAN
Specificity	Human Ribonuclease 7

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

