

Anti-BDNF antibody [3C11] (STJA0002904)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/IF
Host / Source	Mouse
Reactivity	Human/Mouse/Rat/Guinea Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G purified
Dilution Range	ELISA 0.02 to 1 μ g/ml; WB 0.1 to 1 μ g/ml under reducing conditions; IF 2 to 20 μ g/ml. Monoclonal antibody working amount has to be established practically for each particular antigen and assay format.
Formulation	PBS pH 7.4, with 0.1% sodium azide.
Isotype	IgG1
Storage Instruction	Store at -20 to -70 °C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles.

TARGET INFORMATION

Gene ID	627
Gene Symbol	BDNF
UniProt ID	BDNF_HUMAN
Specificity	human, mouse, rat, guinea pig

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

