

Anti-S100-A9 protein antibody [6G10] (STJA0002914)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	NA Protein G purification NA
Dilution Range	ELISA 1-0.2 µg/ml
Formulation	PBS with 50% glycerol
Isotype	IgG1
Storage Instruction	Store at -20 to -70°C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6280
Gene Symbol	S100A9
UniProt ID	S10A9_HUMAN
Specificity	Human S100-A9

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---