

Anti-Zaire Ebolavirus GP protein antibody (STJA0002915)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/IF
Host / Source	Rabbit
Reactivity	EBOV

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Ebolavirus GP protein-affinity purification
Dilution Range	ELISA 0.02 $\mu\text{g/ml}$; WB 0.5 $\mu\text{g/ml}$; IF 0.078 to 5 $\mu\text{g/ml}$ Antibody working titer has to be established practically for each particular antigen and assay format
Formulation	Concentrated ammonium sulphate in PBS pH 7.4
Isotype	IgG
Storage Instruction	Store at +4°C upon receipt. (NH ₄) ₂ SO ₄ precipitate, mix well by pipetting or vortexing prior use.

TARGET INFORMATION

Specificity	Zaire Ebolavirus GP protein
-------------	-----------------------------

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

