

Anti-RNase 8 antibody (STJA0002917)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/IF
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	RNase 8-affinity purified
Dilution Range	ELISA 1 to 2 ng/ml; WB 0, 1 to 0, 5 μ g/ml; IF 1 to 10 μ g/ml Antibody optimal dilution has to be established practically for each antigen and assay format
Formulation	Concentrated ammonium sulphate in PBS pH 7.4
Isotype	IgG
Storage Instruction	Store at +4°C upon receipt. As product is (NH ₄) ₂ SO ₄ precipitate, mix well by pipetting or vortexing prior use

TARGET INFORMATION

Gene ID	122665
Gene Symbol	RNASE8
UniProt ID	RNASE8_HUMAN
Specificity	human RNase 8, cross-reactivity to human RNase 7

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

