

Anti-Recombinant-MANF antibody [3D4] (STJA0002919)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/IF/IHC-C
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	MabSelect affinity chromatography
Dilution Range	ELISA 0.02 to 0.1 μ g/ml; IF 1 to 20 μ g/ml; IHC-C 5 μ g/ml. Conformational antibody, not suitable for Western Blot application. Monoclonal antibody working titer has to be established practically for each particular antigen and assay format
Formulation	PBS pH 7.4, with 0.1% sodium azide
Isotype	IgG2a
Storage Instruction	Store at -20 to -70°C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles

TARGET INFORMATION

Gene ID	7873
Gene Symbol	MANF
UniProt ID	MANF_HUMAN
Specificity	Human MANF

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

