

Anti-Recombinant-SARS-CoV-2 Spike S1 antibody [8E5] (STJA0002920)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA
Host / Source	Human
Reactivity	SARS-CoV-2

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A affinity chromatography following desalting
Dilution Range	Indirect ELISA 0.02 – 0.04 ng/ml
Formulation	PBS pH 7.4
Isotype	IgG1 lambda
Storage Instruction	Store at -20 °C to -70 °C long term. Avoid multiple freeze-thaw cycles.

TARGET INFORMATION

Specificity	SARS-CoV-2 Spike S1
-------------	---------------------

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Coomassie-stained SDS-PAGE analysis of monoclonal antibody to SARS-CoV-2 Spike S1, clone 8E5. 4-12% gradient gel is used for analysis. Lane 1. Monoclonal antibody to SARS-CoV-2 Spike S1, clone 8E5 (-DTT). Lane 2. Size Marker. Lane 3. Monoclonal antibody to SARS-CoV-2 Spike S1, clone 8E5 (+DTT). Lane 4. Size Marker