

Anti-Recombinant-SARS-CoV-2 Spike S2 antibody [1H11] (STJA0002921)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA
Host / Source	Human
Reactivity	SARS-CoV-2

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A affinity chromatography following desalting
Formulation	PBS pH 7.4
Isotype	IgG1 lambda
Storage Instruction	Store at -20 °C to -70 °C long term. Avoid multiple freeze-thaw cycles.

TARGET INFORMATION

Specificity	SARS-CoV-2 Spike S2
-------------	---------------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---