

Anti-Recombinant-HEL Mouse IgG1 antibody (STJA0002923)

GENERAL INFORMATION

Product Type	Primary antibodies
Host / Source	CHO cells
Reactivity	Chicken

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	NA Hitrap Mabselect Sure + Superdex 200 Increase NA
Formulation	PBS pH 7.4
Isotype	mIgG1
Storage Instruction	Store at -20 °C to -70 °C long term. Avoid multiple freeze-thaw cycles.

TARGET INFORMATION

Specificity	Lysozyme from chicken egg white
-------------	---------------------------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------