

Anti-Recombinant-SARS-CoV-2 Spike RBD antibody [71D2] (STJA0002931)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA
Host / Source	CHO cells
Reactivity	SARS-CoV-2

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A affinity chromatography following gel filtration
Dilution Range	Indirect ELISA 0.02 – 0.04 ng/ml
Formulation	PBS pH 7.4
Isotype	IgG kappa
Storage Instruction	The product is stable at 4 °C for several days. Store at -20 °C to -70 °C long term. Avoid multiple freeze-thaw cycles.

TARGET INFORMATION

Specificity	SARS-CoV-2 Spike RBD
-------------	----------------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---