

Anti-Recombinant-SARS-CoV-2 Spike RBD antibody [79C3] (STJA0002933)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA
Host / Source	CHO cells
Reactivity	SARS-CoV-2

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A affinity chromatography following gel filtration
Dilution Range	Indirect ELISA 0.01 – 0.02 ng/ml
Formulation	PBS pH 7.4
Isotype	IgG kappa
Storage Instruction	The product is stable at 4 °C for several days. Store at -20 °C to -70 °C long term. Avoid multiple freeze-thaw cycles.

TARGET INFORMATION

Specificity	SARS-CoV-2 Spike RBD Neutralizing properties: blocks S1 and ACE2 binding in Octet system.
-------------	---

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---