

Anti-Recombinant-interleukin-6 antibody [11D8] (STJA0002937)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA
Host / Source	CHO cells
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	NA Protein A affinity chromatography following desalting NA
Formulation	PBS pH 7.4
Isotype	IgG1
Storage Instruction	Store at -20 °C to -70 °C long term. Avoid multiple freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3569
Gene Symbol	IL6
UniProt ID	IL6_HUMAN
Specificity	Reacts with human IL-6. No cross-reactivity with IL-1 Alpha, IL-1 Beta, IL-2, IL-3, IL-4, IL-8, IFN Gamma and TNF Alpha proteins.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---