

Anti-Recombinant-Procalcitonin antibody [9D8] (STJA0002940)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA
Host / Source	CHO cells
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	NA Protein A affinity chromatography following desalting NA
Formulation	PBS pH 7.4
Isotype	NA IgG1 NA
Storage Instruction	Store at -20 °C to -70 °C long term. Avoid multiple freeze-thaw cycles.

TARGET INFORMATION

Specificity	human PCT
-------------	-----------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---