

## Anti-Recombinant-Procalcitonin antibody [14H11] (STJA0002943)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA
Host / Source	CHO cells
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	NA Protein A affinity chromatography following desalting NA
Formulation	PBS pH 7.4
Isotype	IgG1
Storage Instruction	Store at -20 °C to -70 °C long term. Avoid multiple freeze-thaw cycles.

### TARGET INFORMATION

Specificity	human PCT
-------------	-----------

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---