

Anti-Recombinant-Procalcitonin antibody [14E8] (STJA0002945)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | ELISA |
| Host / Source | CHO cells |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Clonality | Monoclonal |
| Conjugation | Unconjugated |
| Purification | Protein A affinity chromatography following desalting. |
| Formulation | PBS pH 7.4 |
| Isotype | IgG1 |
| Storage Instruction | Store at -20 °C to -70 °C long term. Avoid multiple freeze-thaw cycles. |

TARGET INFORMATION

| | |
|-------------|----------------------------|
| Gene ID | 796 |
| Gene Symbol | CALCA |
| UniProt ID | CALC_HUMAN |
| Specificity | human PCT |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|