

Anti-Recombinant-Cardiac Troponin I antibody [5F8] (STJA0002946)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA
Host / Source	CHO cells
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	MabSelect affinity chromatography following gel filtration
Formulation	50 mM citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ (SPTN)
Isotype	IgG2 lambda
Storage Instruction	Store at -20 °C to -70 °C long term. Avoid multiple freeze-thaw cycles.

TARGET INFORMATION

Specificity	human cardiac troponin I
-------------	--------------------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---