

Anti-Alpha Internexin/NF-66 antibody [1D2] (STJA0003572)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/ICC
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	This antibody was protein g purified culture from supernatant.
Dilution Range	WB 1:2000 IHC 1:1000-1:5000 ICC 1:250-1:5000
Formulation	100 ul in PBS + 10 mM Sodium Azide.
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene Symbol	Ina Uniprot ID= AINX_RAT "
-------------	--

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



IF staining rat cerebellum alpha-internexin in green with anti-calretinin in red STJA0003572



IHC staining rat CNS cells with Anti-alpha-internexin neuronal processes in red STJA0003572



Western blot 66 kDa alpha internexin protein in rat cortex lysate STJA0003572