

Anti-CALB1 antibody (Full-Length) (STJA0003587)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC
Host / Source	Mouse
Reactivity	Bovine/Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	This antibody was protein g purified culture from supernatant.
Dilution Range	WB 1:2000 IHC 1:1000-1:5000
Formulation	100 ul in PBS + 50% Glycerol and 5 mM Sodium Azide
Isotype	IgG2a
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	793
Gene Symbol	CALB1
UniProt ID	CALB1_HUMAN
Immunogen Region	Full-Length

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



IHC staining rat cerebellum using anti-calbindin antibody Purkinje cell dendrites STJA0003587



Western blot immunolabeling STJA0003587
 28 kDa calbindin protein
 in rat cerebellar lysate



IHC staining with calbindin antibody in rat cerebellum Purkinje cell dendrites in red STJA0003587