

Anti-DOPA Decarboxylase antibody (STJA0003612)

GENERAL INFORMATION


Product Type	Primary antibodies
Applications	WB
Host / Source	Rabbit
Reactivity	Bovine/Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	This antibody was antigen affinity purified from pooled serum.
Dilution Range	WB 1:1000
Formulation	100 μ l in 10 mM HEPES (pH 7.5) , 150 mM NaCl, 100 μ g per ml BSA and 50% Glycerol.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



IHC staining mouse brain section DDC immunolabeling using antibody STJA0003612



Western blot 55 kDa DDC protein immunolabeling in rat adrenal medulla STJA0003612