

Anti-Dopamine Beta-Hydroxylase antibody (C-Term) (STJA0003614)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Sheep
Reactivity	Human/Mouse/Non-Human Primates

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	This antibody was antigen affinity purified from pooled serum.
Dilution Range	WB 1:1000
Formulation	100 μ l in 10 mM HEPES (pH 7.5) , 150 mM NaCl, 100 μ g per ml BSA and 50% Glycerol.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1621
Gene Symbol	DBH
UniProt ID	DOPO_HUMAN
Immunogen Region	C-Term

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



antibody Western blot 75 kDa DBH protein in human adrenal medulla lysate STJA0003614