

## Anti-GRIA2/3 antibody (C-Term) (STJA0003664)

### GENERAL INFORMATION

|               |                                   |
|---------------|-----------------------------------|
| Product Type  | Primary antibodies                |
| Applications  | WB                                |
| Host / Source | Rabbit                            |
| Reactivity    | Rat/Chicken/Human/Mouse/Zebrafish |

### PRODUCT PROPERTIES


|                     |                                                                                                |
|---------------------|------------------------------------------------------------------------------------------------|
| Clonality           | Polyclonal                                                                                     |
| Conjugation         | Unconjugated                                                                                   |
| Purification        | This antibody was antigen affinity purified from pooled serum.                                 |
| Dilution Range      | WB 1:1000                                                                                      |
| Formulation         | 100 $\mu$ l in 10 mM HEPES (pH 7.5) , 150 mM NaCl, 100 $\mu$ g per ml BSA and 50% Glycerol.    |
| Isotype             | IgG                                                                                            |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

### TARGET INFORMATION

|                  |                           |
|------------------|---------------------------|
| Gene ID          | <a href="#">29627</a>     |
| Gene Symbol      | <a href="#">Gria2</a>     |
| UniProt ID       | <a href="#">GRIA2_RAT</a> |
| Immunogen Region | C-Term                    |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|



Western blot GluR2-3 protein 100 kDa in rat hippocampal lysate STJA0003664