

Anti-LMNA antibody (Full-Length) [4C4] (STJA0003674)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/ICC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	This antibody was protein g purified culture from supernatant.
Dilution Range	WB 1:5000 IHC 1:1000 ICC 1:1000
Formulation	100 ul in PBS + 10 mM Sodium Azide.
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	4000
Gene Symbol	LMNA
UniProt ID	LMNA_HUMAN
Immunogen Region	Full-Length

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC staining HeLa cells with anti-Lamin A-C antibody nuclear lamina localization STJA0003674



Western blot lamin A-C protein bands at 64 kDa and 74 kDa in HeLa cell lysate STJA0003674



ICC staining HeLa cells with anti-Lamin A-C red nuclear lamina localization STJA0003674