

Anti-LAMP1 antibody (Full-Length) [5H6] (STJA0003675)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/ICC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	This antibody was protein g purified culture from supernatant.
Dilution Range	WB 1:5000 IHC 1:2000 ICC 1:500-1:2000
Formulation	100 ul in PBS + 10 mM Sodium Azide.
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	3916
Gene Symbol	LAMP1
UniProt ID	LAMP1_HUMAN
Immunogen Region	Full-Length

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC staining HeLa cells with anti-LAMP1 antibody
lysosomal vesicular pattern in red STJA0003675



Western blot 100 kDa LAMP1 protein detection in
HeLa cell lysate STJA0003675



IHC staining HeLa cells with anti-LAMP1 antibody
lysosomal cytoplasmic labeling STJA0003675