

Anti-Phospho-STXBP2-Ser 513 antibody (STJA0003694)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC
Host / Source	Rabbit
Reactivity	Rat/Bovine/Canine/Chicken/Human/Mouse/Non-Human Primates/Xenopus

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	This antibody was antigen affinity purified from pooled serum.
Dilution Range	WB 1:1000 ICC 1:400
Formulation	100 μ l in 10 mM HEPES (pH 7.5) , 150 mM NaCl, 100 μ g per ml BSA and 50% Glycerol.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	81804
Gene Symbol	Stxbp2
UniProt ID	STXB2_RAT

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot immunolabeling of Munc-18 in rat cortical lysate targeting phospho-Munc-18 Ser513 in red. STJA0003694



IHC staining rat cortical neurons with antibody targeting phospho-Munc-18 Ser513 in red. STJA0003694