

Anti-NES antibody [4D11] (STJA0003696)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/ICC
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	This antibody was protein g purified culture from supernatant.
Dilution Range	WB 1:2000 IHC 1:500 ICC 1:500
Formulation	100 ul in PBS + 50% Glycerol and 5 mm Sodium Azide
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	10763
Gene Symbol	NES
UniProt ID	NEST_HUMAN

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Western blot nestin doublet immunolabeling in neonatal rat brain lysate STJA0003696



ICC staining with anti-nestin antibody rat cortical neurons and glia with MAP2 co-stain STJA0003696



ICC staining rat neurons and glia with anti-nestin antibody red and anti-vimentin green STJA0003696



IHC staining embryonic mouse cortex with Nestin antibody green and DAPI nuclear stain blue STJA0003696