

Anti-NEFM antibody (C-Term) [3H11] (STJA0003701)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/ICC
Host / Source	Mouse
Reactivity	Bovine/Chicken/Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	This antibody was protein g purified culture from supernatant.
Dilution Range	WB 1:1000 IHC 1:500-1:5000 ICC 1:100-1:500
Formulation	100 ul in PBS + 10 mM Sodium Azide.
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	24588
Gene Symbol	Nefm
UniProt ID	NFM_RAT
Immunogen Region	C-Term

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot immunolabeling STJA0003701 145 kDa NF-M rat cortical



protein anti-Neurofilament-M IHC staining rat frontal cortex pyramidal neuron cell bodies



ICC staining cultured rat neurons with antibody NF-M labeling in green STJA0003701