

Anti-GRIN2A antibody (C-Term) (STJA0003711)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/IP
Host / Source	Rabbit
Reactivity	Human/Mouse/Rabbit/Rat

PRODUCT PROPERTIES



Clonality	Polyclonal
Conjugation	Unconjugated
Purification	This antibody was antigen affinity purified from pooled serum.
Dilution Range	WB 1:1000 IHC 1:1000 IP 3 µl per 200 µg lysate
Formulation	Lyophilized
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	24409
Gene Symbol	Grin2a
UniProt ID	NMDE1_RAT
Immunogen Region	C-Term

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot immunolabeling STJA0003711 NR2A subunit 180 kDa in rat hippocampal lysate IHC staining with NR2A in rabbit retina photoreceptors and plexiform layers STJA0003711